

Empire Diagnostic Solutions, Inc.

1090 Coney Island Ave 1st Floor Brooklyn, NY 11230
Tel: 347-715-4665, Fax: 347-715-4668
www.edslive.com



Important information you should know during your 72 hour Ambulatory EEG Test (AEEG)

- Complete the “Event Calendar”.
- **NO** showering while the test is in progress.
 - NOTE: You may wash from the neck down using a wash cloth. Be careful not to get any water onto the machine.
- **NO** gum chewing during the test.
- **NO** tampering of any kind with the machine.
- **DO NOT** drop the machine at any time during the test.
- Once test is complete the electrodes may be removed as follows:
 - Acetone or nail polish remover may be used by soaking each electrode with a piece of gauze or cotton for 30 seconds. Once electrodes have been removed you may wash with warm water and shampoo.
 - WARNING: Make sure acetone or nail polish remover does not get into your eyes. If that happens then flush with cold water.
 - Please, remove the electrodes on _____.
 - Return the equipment, electrodes and completed event calendar on the same day between 9 AM to 5 PM at _____.

In case of emergency or questions during the 72 hour test, please call our service at 800-373-1843 and request to speak to the manager.

IMPORTANT NOTE: The patient must remember to return the backpack with all equipment to your office. We will pick up the backpack within 24 hours.